


TASK Lingual Orthodontic Instruments

LINGUAL BOND REMOVING PLIERS (for posterior teeth)

- Adhesive material at hard-to-reach area for conventional remover can be easily removed.
- Tip for replacement is available.
- Pliers are autoclavable.



LINGUAL BOND REMOVING PLIERS (for posterior teeth)

Model # : P-700

Price : **\$141.30**

EXCHANGEABLE PARTS

METAL TIP

Model # : P-701

Price : **\$50.60**

PEEK TIP (3 / bag)

Model # : 9836

Price : **\$13.00**



17972 Sky Park Circle, Suite J, Irvine, CA 92614
Phone: 949-474-0176 Fax: 949-474-1736
E-mail: charlie@headdental.com
website: www.headdental.com