

FINENESS

Ultra Thin, Light Weight Resin Forming Instrument "FINENESS" Lines

The instrument coated by titanium nitride, is comfortable to use. Ultra thin tip achieves easy access to fill and form resin in area that has traditionally been difficult to reach.

Ultra Thin
0.17mm

Light
Weight
15g

Hollow stainless steel handle is adopted. You have great selection among all the lines of easy-to-use instruments.

Ultra Thin, Light Weight 4 Different Fineness

Fineness 1



Fineness 2



Fineness 3



Fineness MI



Fineness

Please see the back side for product specification of Resin Forming Instrument "FINENESS" Lines.