



Laboratory Micromotor Set

LMS-20S

■ Compact Design with Light Weight

The total weight, less than 15 oz., allows you to carry the unit anywhere you like. It is the perfect tool to adjust / polish denture outside office.

■ Simple Design

The control panel is designed simple for its easy use.

Suggested
Retail Price
of LMS-20S :
\$458

■ LMS-20S includes

LBM-100	E-type Micromotor
CNT-001	Controller
AC Adaptor	
ES-30A	E-type Straight Handpiece

■ Specification

• Speed :	2,000 - 20,000 rpm	• Weight :	LBM-100 5.5 oz
• Torque :	0.22 lbf·in		CNT-001 1.6 oz
• Electric Power Source :	AC100 - 240V, 50/60Hz		AC Adaptor 5.5 oz
			ES-30A 1.8 oz



The micromotor can be powered by portable battery, besides AC adaptor. The use of portable battery makes your mobile area boundless.

* Output voltage of portable battery can affect the maximum speed of micromotor.

Battery Recommended

- Output Voltage :
16 - 24V
- Output Current :
1.5 - 3.0A



Distributed by