

AUTO STRIPPING BUR - Double Sided



Double Sided Strip is a thin strip with diamond coating on two sides.

DOUBLE SIDED

COLOR	PRODUCT CODE	DENSITY OF DIAMOND COATING ON THE BLADE (Unit : μm)	THICKNESS OF BLADE (Unit : mm)
P	001-015H	15 (FINE)	0.08
P	001-025H	25	0.12
P	001-040H	40	0.15
P	001-060H	60	0.19
P	001-090H	90 (COARSE)	0.29

1 μm + 1/1000 mm.
Thickness of Blade is not same size as a gap between the teeth to strip.